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Stand: Juli 2015

## Processing instructions

in accordance with DIN EN ISO 17664:2004 (re-usable medical devices)

### Product(s):

All re-usable surgical instruments supplied by Albert Heiss GmbH & Co. KG which can be defined as fixed constructions (no moving parts), as assemblies (partially movable surgical instruments which can be disassembled) or as joined constructions (simple jointed constructions which cannot be disassembled).

### WARNING NOTES:

Special care is necessary when cleaning narrow gaps between parts, grooves, through-holes and blind holes, hinge and screwed constructions in addition to parts lying on or moving over each other.

### Limitation of reprocessing:

Frequent treatment has little or no effect on these instruments. Therefore the end of the product's useful life depends mainly on mechanical treatment and stress.

### Area of use

Ophthalmic surgical instruments are used in ophthalmic surgery. They are specifically suitable for restricted spaces. They are used for cutting, grasping, dissecting, dilating, draining, aspirating, etc. of vessels and tissue.

### Models and sizes supplied

The current product range is described in the Albert Heiss catalogue.

### General notes

These instructions for processing are applicable for all products manufactured by Albert Heiss GmbH & Co. KG. Supplementary instructions, e.g. on disassembly or the use of flushing adaptors are enclosed with each product and should be seen as supplementary to these instructions.

General deviations from these instructions regarding other sterilization procedures or during manual or mechanical cleaning and disinfection are the user's responsibility. Albert Heiss GmbH & Co. KG accepts no liability for product damage arising due to deviations from these instructions. Our General Terms and Conditions apply in addition. Albert Heiss GmbH & Co. KG reserves the right to make changes to these instructions on the basis of more recent knowledge.

The most up-to-date version of these instructions is available at [www.albert-heiss.de](http://www.albert-heiss.de) and can also be ordered from Albert Heiss GmbH & Co. KG.

### Basic principles

All instruments are supplied in a non-sterile state and must be cleaned/disinfected and sterilized before use. Efficient cleaning/disinfection is the requirement for effective sterilization of the product.

The legal provisions applicable in each country and the hygiene regulations of the hospital or clinic must be observed in addition to these instructions.

Because of the intended use (contact with blood), we recommend categorizing instruments with complex geometries (e.g. lumen and hinge instruments) as "critical B" in accordance with RKI guidelines. Instruments with simple geometries (e.g. knives, lances, forceps) may be categorized either as "critical A" or "semi-critical A" as long as there is no contact with blood.

## PROCESSING INSTRUCTIONS

### A: Storage and transport:

It is recommended that the processing of an instrument is carried out immediately after its use, in particular after contact with foreign matter.

### B: Preparation for cleaning:

For effective reprocessing, pre-treatment should take place no later than 30 minutes after the end of the operation.

Ensure unobstructed flow in lumen instruments (e.g. cannulas, handpieces). Cleaning/disinfection should be started within the next 2 hours.

The following pre-treatment steps must be performed regardless of the type of cleaning/disinfection process (manual or mechanical) subsequently used:

### Pre-treatment sequence:

1. Disconnect all detachable connections.
2. Remove all visible contaminants with a moist cloth.
3. Flush any connectors for aspiration and/or irrigation lines at least 5 times with 10 ml deionized or distilled water ensuring that the water flows in one direction only.
4. Aldehyde-free disinfectant may be added to the flushing water. Subsequently, the lines must be flushed an additional 5 times with 10 ml of deionized or distilled water.

### C: Ultrasonic cleaning

Heavily soiled instruments may additionally be pre-cleaned in an ultrasonic bath. In this case, care must be taken to adhere to the exposure time and concentrations specified by the manufacturer of the cleaning solution.

### General notes:

- The cleaning solution must completely cover the products.
- Instruments with hinges or joints, e.g. scissors, must be treated in an open state
- Only perforated trays may be used to store the instruments in order not to impair the effect of the ultrasound.
- Replace the cleaning solution regularly depending on the conditions of use.

### D: Mechanical cleaning:

#### General notes on the cleaning and disinfection equipment (CDE) to be used:

- The CDE must meet the requirements of DIN EN ISO 15883/AAMI ST15883 and have a verified effectiveness.
  - Flushing adaptors are required to allow the connection of lumen instruments (e.g. cannulas, handpieces). These adaptors must have flushing pressures which have been validated user-specifically as being sufficient and reproducible.
  - Thermal disinfection is preferable to chemical disinfection. When disinfecting with chemicals there is a risk of disinfectant residues remaining on the products if final rinsing is inadequate.
  - The program used must be appropriate for the products and must include a sufficient number of rinsing cycles.
  - Only sterile or low-bacteria water must be used.
  - The air used for drying the instruments must be filtered. The air must also meet the hospital's requirements for compressed air.
- The following must be taken into account when selecting cleaning agents:
- The cleaning agent used must be fundamentally suitable for cleaning of the products.
  - A suitable disinfectant with verified effectiveness should be used if thermal disinfection is not carried out. The disinfectant must be compatible with the detergent used.
  - If an alkaline detergent is used, adequate neutralization in accordance with the manufacturer's instructions must be performed.
  - It is imperative to adhere to the detergent and disinfectant concentrations specified by the manufacturer.

### Mechanical cleaning sequence:

1. Place the instruments in a perforated sterilizing tray and fix the products in place so that the effectiveness of cleaning/disinfection is not impaired; ensure the products do not touch each other during the cleaning/disinfection process.
2. Place the perforated sterilizing tray containing the products in the CDE.
3. Attention: The manufacturer's instructions must be followed when stacking several perforated sterilizing trays one on top of the other.
4. Connect any lumen instruments (e.g., cannulas, handpieces) to the flushing adaptors of the CDE and, if possible, close any remaining unused flushing adaptors.
5. Attention: Always check instruments with a lumen for unobstructed flow before the cleaning/disinfection process otherwise the effectiveness of the cleaning process may be called into question.
6. Start the program on the CDE.
7. When the program is finished, take the perforated sterilizing tray holding the products out of the CDE.
8. It may be necessary to dry instruments with a lumen using filtered compressed air.
9. The products must be packaged immediately after being taken out of the CDE.

### E: Manual cleaning:

Due to its significantly higher effectiveness, a mechanical process is always preferable to the manual process for cleaning/disinfection.

### Procedure:

1. Rinse the instrument to remove surface contaminants
2. Place the instrument in a deep-cleaning disinfectant so that all surfaces, cavities, lumina and openings are covered. Observe the disinfectant manufacturer's instructions.
3. Instruments with hinges or joints must be cleaned in both the open and also the closed position.
4. Remove contaminants adhering to the outside using a soft synthetic brush. Never use abrasive cleaning agents or metal brushes.
5. Thoroughly rinse all instruments with a lumen (e.g. cannulas, handpieces) with cleaning solution 5 - 10 times at the beginning and end of the exposure time.
6. Carry out final rinsing with distilled or demineralized water. Observe the cleaning agent manufacturer's instructions
7. Place the instrument in a disinfectant so that all surfaces, cavities, lumina and openings are covered. Observe the disinfectant manufacturer's instructions.
8. Instruments with hinges or joints must be disinfected in both the open and also the closed position.
9. Thoroughly rinse all instruments with a lumen (e.g. cannulas, handpieces) with cleaning agent 5 - 10 times at the beginning and end of the exposure time.
10. Carry out final rinsing at least 10 times with distilled or demineralized water. Observe the disinfectant manufacturer's instructions.
11. Repeat the complete cleaning/disinfection process if the last rinsing solution is not clear or if contaminants are still visible on the product.
12. Subsequently dry the instrument with an absorbent, soft, lint-free cloth.
13. Dry the products, including all lumina and canals, using filtered compressed air.

### F: Inspection and functional test:

- Allow instrument to cool to ambient temperature so that there is no occurrence of metal abrasion which can lead to corrosion.
- Lightly oil movable parts (e.g. hinges and locks) with sterilizable, steam-permeable oil (LubriPEN, Art Nr. H-7123).
- Check instrument after every cleaning and disinfection procedure for: cleanliness, function and damage,
- sort out and replace any damaged and defective instrument.

### G: Preparation for sterilization:

It is necessary to ensure that hinges and joints are open and that water can drain out of narrow gaps, grooves, through-holes and blind holes as well as hinged, jointed and screwed constructions. Sterilization with the instruments uncovered or in minimum contact is recommended.

### H: Recommended sterilization:

The times and temperatures shown below are minimum requirements and must be complied with. If these values have to be lowered, the user must revalidate them.

In principle, the times and temperatures may be exceeded. However, this leads to increased stress on the material and causes the products to age faster.

Only products which have been cleaned and disinfected may be sterilized.

The sterilization processes shown below must be used for sterilization. The user must revalidate the process if there is any deviation from this requirement.

- Fractionated vacuum process with adequate drying of the product
- Steam sterilizer in accordance with DIN EN 13060 or DIN EN 285 and validated in accordance with DIN EN (ANSI AAMI) ISO 17665-1
- Maximum sterilization temperature 137 °C (plus tolerance range in accordance with DIN EN (ANSI AAMI) ISO 17665-1).
- Sterilization times:
  - at least 20 minutes at 121 °C, alternatively
  - at least 5 minutes at 132/134 °C<sup>\*</sup>

<sup>\*</sup>) The Robert Koch Institute recommends a sterilization time of at least 5 minutes at 134 °C.

This recommendation requires the use of an alkaline cleaning agent with a pH-value which is greater than 10. Should an alkaline cleaning agent not be used, then the sterilization time must be increased to 18 minutes.

### I: Packaging:

1. Store instruments with a fine point and/or microsurgery instruments in appropriate storage aids

2. These must be placed in a suitable sterilization container for sterilization of the products. This, together with the filter material, must meet the following criteria:

- be in accordance with the standards DIN EN (ANSI AAMI) ISO 11607 (DIN EN 868-8)
- be suitable for steam sterilization (temperature-resistant up to 141 °C, adequate steam permeability)
- be serviced regularly

If disposable sterilization packages are used, these must also be in accordance with DIN EN (ANSI AAMI) ISO 11607 and suitable for steam sterilization (temperature-resistant up to 141 °C, adequate steam permeability).

### Storage:

Sterilized products should be stored in a dry environment. No other special conditions are necessary as far as storage is concerned.

### Supplementary information:

Do not exceed the sterilizer's maximum load when sterilizing several instruments in one sterilization cycle.

- Products must be stored appropriately during sterilization (e.g. not placed on top of each other and fixed in place with sterilization strips or sterilization plates).
- Products should be sterilized in an "open" state
- Only products suitable for sterilization are to be used in the care of hinges and joints of instruments.
- Titanium instruments may become discolored when using hydrogen peroxide H<sub>2</sub>O<sub>2</sub>. These discolorations can be attributed to changes in the thickness of the oxide layer and do not affect the quality of the instruments.

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